

## Jeio Tech Water Baths (BW-B Series)

**Brand:** Jeio Tech

Constructed with a corrosion-resistant stainless steel tank, Jeio Tech's Water Baths (BW-B Series) are an economical range of water baths that provide stable temperature control for samples.

Each unit comes with an easy and fast sloped control panel, a safety alarm setting, and splash-proof keypad.

### Features:

#### Performance

- Temperature range from ambient +7°C to 100°C.
- On/off temperature control.

#### Convenience

- Simple and quick control panel.
- Corrosion resistant stainless steel construction allowing other bath fluids as well as water and silicone fluids.
- Low-profile design with sloped control panel for easier access and better safety.
- Variety of accessories are available to meet application demands.

#### Safety

- Isolated design of heater and sensor by a baffle plate to protect users and samples against direct contact.
- Complete safety protection system with warning alarms. Over-temperature, over-current protection.
- Splash-proof keypad.



## Specifications:

Model(s)		BW-05B	BW-10B	BW-20B
Product Code(s)		AAH46012K	AAH46022K	AAH46032K
Bath Volume (litres)		3.5	11.5	20
Working Temperature Range (°C)		Amb. +7 to 100		
Temperature Stability (°C)		± 0.5		
Dimensions	Bath Opening/Depth (w x l, d)	240 x 136, 115mm	300 x 240, 165mm	498 x 300, 165mm
	External (w x d x h)	303 x 242 x 245mm	357 x 330 x 290mm	564 x 392 x 290mm
	Net Weight (kg)	6.5	10	18.5
Electrical Requirements (230V)		50Hz, 3A	50Hz, 4.3A	50Hz, 8.7A

Unit 3, Tower Business Park  
 Warpsgrove Lane  
 Chalgrove  
 Oxfordshire  
 OX44 7XZ

Tel: (+44) 01865 400321

Email: [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

Website: [www.medlinescientific.com](http://www.medlinescientific.com)