

Medline Water Baths (MDBS Series)

Brand: Medline Scientific

Manufactured from polycarbonate with a safety thermostat, Medline's Water Baths (MDBS Series) provide optimal temperature control for your samples.



Specifications:

Product Code(s)		MDBS-3	MDBS-7	MDBS-14	MDBS-30
Bath Volume (litres)		3	7	14	30
Temperature	Range (°C)	Amb. +5 to 100			
	Accuracy (°C)	± 0.3 (at 37°C)			
Dimensions	Bath Opening/ Depth (w x l, d)	100 x 170, 100mm	160 x 300, 110mm	210 x 350, 160mm	320 x 510, 180mm
	External (w x d x h)	150 x 310 x 240mm	220 x 440 x 240mm	280 x 500 x 300mm	380 x 660 x 240mm
Weight (kg)		2	3.5	6.5	11.2

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com