

 CAPP APS  
LANDBRUGSVEJ 10  
DK-5260 ODENSE  
DENMARK

TEL. (+45) 6613 6140  
FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK  
WEB: WWW.CAPP.DK

**CAPP**  
D E N M A R K

## Capp *Rondo*

M A G N E T I C S T I R R E R

The Capp Rondo Magnetic Stirrer is a compact device that works at a maximum speed of 1500 rpm. It is manufactured with robust stainless steel top surface to resist moderate contamination of harsh chemicals. The adjustable speed and timer function can be set and monitored on a modern real-time digital display.

### Not just a simple Stirrer:

Adjust speed and timer.

Last run memory feature.

Pulse mode works for 30 seconds on each direction, for an improved mixing.

Water spill resistant – IP 65 compliant.

### Innovative design:

Slim profile – only 16.5mm thick.

Touch control pad and digital display.

Flat surface for ease of cleaning.

Maintenance free.

**Important note:** to maximize the performance of your Capp Rondo Magnetic Stirrer, Capp recommends that you use our original Teflon coated, 25mm high strength stirring bar. One unit is supplied with each instrument. Contact your local distributor to order more units.

Specifications	Capp Rondo Magnetic Stirrer
Maximum speed	1500 rpm
Speed setting	15 to 1500 rpm (increments of 50 rpm)
Load capacity	800 mL
Recommended stir bar length	25 mm
Dimensions	124 x 191 x 16.5mm
Weight	600 grams
Pulse function	YES
Timer function	1 to 99 minutes and infinite mode



## Ordering Info

Cat No. Description

CRS-15X Capp Rondo Magnetic Stirrer 1500rpm w/ 1 stirring bar