

CAPP *Tronic*

ELECTRONIC PIPETTES



The **CAPP Tronic** is certainly one of the easiest electronic pipettes to set and use.

Having only a one-side programming manual, familiarizing yourself with **CAPP Tronic's** functions takes just a matter of minutes.

Apart from standard pipetting, you can choose between the most widely used pipetting modes - including reverse pipetting, multi dispensing, mixing and sequential dispensing.

Up to nine individual protocols can be stored thanks to the extensive memory capacity.

Supplied with one long-life and easily exchangeable lithium battery, the **CAPP Tronic** allows for uninterrupted pipetting without the need to stop work in order to recharge.

The O-rings on the tip base provides an excellent fit with a great variety of universal fitting tips.



Ordering Information

Cat. No.	Description	Inaccuracy %	Imprecision %	Colour
E20-8BA	CAPPTronic electronic pipette, 8-channel, 2-20 μ L	1.00	0.30	●
E200-8BA	CAPPTronic electronic pipette, 8-channel, 10-200 μ L	0.60	0.15	●
E1200-8BA	CAPPTronic electronic pipette, 8-channel, 100-1200 μ L	0.50	0.15	●